

SOFIA™ 88

Neurovascular Support Catheter



 **TERUMO**
NEURO



SOFIA™ 88 Neurovascular Support Catheter

Product Code	Product Name	Catheter Size (F)	ID (F / in / mm)	Distal OD (F / in / mm)	Proximal OD (F / in / mm)	Total Working Length (cm)	Distal/Proximal Length (cm)
DA8115ST	SOFIA™ 88	8	6.7 / 0.088 / 2.24	7.8 / 0.102 / 2.6	8.2 / 0.108 / 2.7	115	14 / 101

Compatible with **SOFIA Flow Plus 6F 131cm** Aspiration Catheter.*
Compatible with 8F Groin/Guiding Sheaths.*



SOFIA™ Optimized Braid and Coil Construction

Indications for Use

The SOFIA 88 Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA 88 Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA 88 Catheter is not intended for use in coronary arteries.

For Healthcare Professional Intended Use Only.

Please refer to IFU for the full list of risks, contraindications, warnings, and precautions.

RX Only: Federal (United States) law restricts this device to sale by or on the order of a physician.

*Data on file at Terumo Neuro.

SOFIA™ is a trademark of MicroVention, Inc., registered in the United States and other jurisdictions. All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Refer to Instructions for Use, contraindications and warnings for additional information.



**Terumo Neuro Worldwide
Innovation Center**
(MicroVention, Inc.)

35 Enterprise
Aliso Viejo, CA 92656
USA

PH +1.714.247.8000
PH 1.800.990.8368